

Product datasheet for **TB424294**

UBE1C (UBA3) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of UBE1C (NM_198197), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224294
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	UBE1C (UBA3) Mouse Monoclonal Antibody [Clone ID: OTI7H6] (TA811903)
ACCN:	NM_198197 , NP_937840
Synonyms:	UBA3, hUba3, MGC22384, DKFZp566j164
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_937840
Locus ID:	9039
Cytogenetics:	3p14.1
Protein Pathways:	Ubiquitin mediated proteolysis



[View online »](#)