

## Product datasheet for **TB424349**

### GRK4 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GRK4 (NM_001004056), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224349
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	GRK4 Mouse Monoclonal Antibody [Clone ID: OT19A3] (TA806069)
ACCN:	<a href="#">NM_001004056</a> , <a href="#">NP_001004056</a>
Synonyms:	GPRK2L; GPRK4; GRK4a; IT11
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001004056</a>
Locus ID:	2868
Cytogenetics:	4p16.3
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Chemokine signaling pathway, Endocytosis



[View online »](#)