

## Product datasheet for **TB424414**

### C20orf30 (TMEM230) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TMEM230 (NM_001009925), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224414
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	C20orf30 (TMEM230) Mouse Monoclonal Antibody [Clone ID: OTI2C8] (TA504888)
ACCN:	<a href="#">NM_001009925</a> , <a href="#">NP_001009925</a>
Synonyms:	C20orf30; dj1116H23.2.1; HSPC274
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001009925</a>
Locus ID:	29058
Cytogenetics:	20p13-p12.3
Protein Families:	Transmembrane



[View online »](#)