

## Product datasheet for **TB424760**

### SETD2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SETD2 (NM_014159) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224760
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	SETD2 Mouse Monoclonal Antibody [Clone ID: OTI1E1] (TA804377)
ACCN:	<a href="#">NM_014159</a> , <a href="#">NP_054878</a>
Synonyms:	HBP231; HIF-1; HIP-1; HSPC069; HYPB; KMT3A; LLS; p231HBP; SET2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_054878</a>
Locus ID:	29072
Cytogenetics:	3p21.31
Protein Families:	Druggable Genome
Protein Pathways:	Lysine degradation



[View online »](#)