

## Product datasheet for **TB424888**

### H3FL (HIST1H3B) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HIST1H3B (NM_003537) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC224888
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<a href="#">NM_003537</a> , <a href="#">NP_003528</a>
Synonyms:	H3I; H3C1; H3C3; H3C4; H3C6; H3C7; H3C8; H3C10; H3C11; H3C12; H3FL; HIST1H3B
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_003528</a>
Locus ID:	8358
Cytogenetics:	6p22.2
Protein Pathways:	Systemic lupus erythematosus



[View online »](#)