

Product datasheet for **TB425118**

DPP1 (CTSC) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CTSC (NM_001114173), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC225118
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001114173 , NP_001107645
Synonyms:	CPPI; DPP-I; DPP1; DPPI; HMS; JP; JPD; PALS; PDON1; PLS
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001107645
Locus ID:	1075
Cytogenetics:	11q14.2
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Lysosome



[View online »](#)