

## Product datasheet for **TB425317**

### EIF4E CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of EIF4E (NM_001130679), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC225317
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	EIF4E Mouse Monoclonal Antibody [Clone ID: OTI5D11] (TA503042)
ACCN:	<u><a href="#">NM_001130679</a></u> , <u><a href="#">NP_001124151</a></u>
Synonyms:	AUTS19; CBP; eIF-4E; EIF4E1; EIF4EL1; EIF4F
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001124151</a></u>
Locus ID:	1977
Cytogenetics:	4q23
Protein Pathways:	Insulin signaling pathway, mTOR signaling pathway



[View online »](#)