

## Product datasheet for **TB425700**

### Cbl c (CBLC) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CBLC (NM_001130852), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC225700
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Cbl c (CBLC) Mouse Monoclonal Antibody [Clone ID: OTI3D9] (TA505052)
ACCN:	<a href="#">NM_001130852</a> , <a href="#">NP_001124324</a>
Synonyms:	CBL-3; CBL-SL; RNF57
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001124324</a>
Locus ID:	23624
Cytogenetics:	19q13.32
Protein Families:	Druggable Genome
Protein Pathways:	Chronic myeloid leukemia, Endocytosis, ErbB signaling pathway, Insulin signaling pathway, Jak-STAT signaling pathway, Pathways in cancer, T cell receptor signaling pathway, Ubiquitin mediated proteolysis



[View online »](#)