

## Product datasheet for **TB425857**

### MECP2 CytoSection

#### Product data:

|                                       |                                                                                                                                                                                                                         |
|---------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:                         | CytoSections                                                                                                                                                                                                            |
| Description:                          | Transient overexpression of MECP2 (NM_001110792), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining                                                                               |
| Species:                              | Human                                                                                                                                                                                                                   |
| Expression Host:                      | HEK293T                                                                                                                                                                                                                 |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC225857                                                                                                                                                                                                  |
| Tag:                                  | C-MYC/DDK                                                                                                                                                                                                               |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)                                                                                                                                      |
| Target Detection Antibodies:          | MECP2 Mouse Monoclonal Antibody [Clone ID: OT12F1] (TA812397)                                                                                                                                                           |
| ACCN:                                 | <a href="#">NM_001110792</a> , <a href="#">NP_001104262</a>                                                                                                                                                             |
| Synonyms:                             | AUTSX3; MRX16; MRX79; MRXS13; MRXSL; PPMX; RS; RTS; RTT                                                                                                                                                                 |
| Storage:                              | Room Temperature, or 2-8°C for long term storage                                                                                                                                                                        |
| Stability:                            | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.                                                                                                                 |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.                                                                                                                   |
| RefSeq:                               | <a href="#">NP_001104262</a>                                                                                                                                                                                            |
| Locus ID:                             | 4204                                                                                                                                                                                                                    |
| Cytogenetics:                         | Xq28                                                                                                                                                                                                                    |
| Protein Families:                     | Druggable Genome                                                                                                                                                                                                        |



[View online »](#)