

## Product datasheet for **TB425863**

### ZBTB37 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ZBTB37 (NM_001122770), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC225863
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ZBTB37 Mouse Monoclonal Antibody [Clone ID: OTI1D7] (TA803616)
ACCN:	<u>NM_001122770</u> , <u>NP_001116242</u>
Synonyms:	D430004I08Rik; ZNF908
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001116242</u>
Locus ID:	84614
Cytogenetics:	1q25.1
Protein Families:	Transcription Factors



[View online »](#)