

Product datasheet for **TB425918**

ASH2L CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ASH2L (NM_001105214), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC225918
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ASH2L Mouse Monoclonal Antibody [Clone ID: OTI2H3] (TA810988)
ACCN:	NM_001105214 , NP_001098684
Synonyms:	ASH2; ASH2L1; ASH2L2; Bre2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001098684
Locus ID:	9070
Cytogenetics:	8p11.23
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)