

Product datasheet for **TB426113**

PSD95 (DLG4) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DLG4 (NM_001128827), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC226113
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PSD95 (DLG4) Mouse Monoclonal Antibody [Clone ID: OT112A5] (TA809015)
ACCN:	NM_001128827 , NP_001122299
Synonyms:	MRD62; PSD95; SAP-90; SAP90
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001122299
Locus ID:	1742
Cytogenetics:	17p13.1
Protein Families:	Druggable Genome
Protein Pathways:	Huntington's disease



[View online »](#)