

## Product datasheet for **TB426433**

### Laminin alpha 4 (LAMA4) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of LAMA4 (NM_001105207), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC226433
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Laminin alpha 4 (LAMA4) Mouse Monoclonal Antibody [Clone ID: OTI5G8] (TA805721)
ACCN:	<a href="#">NM_001105207</a> , <a href="#">NP_001098677</a>
Synonyms:	CMD1JJ; LAMA3; LAMA4*-1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001098677</a>
Locus ID:	3910
Cytogenetics:	6q21
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	ECM-receptor interaction, Focal adhesion, Pathways in cancer, Small cell lung cancer



[View online »](#)