

Product datasheet for **TB427419**

Troponin I fast skeletal muscle (TNNI2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TNNI2 (NM_001145829), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC227419
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Troponin I fast skeletal muscle (TNNI2) Mouse Monoclonal Antibody [Clone ID: OT13A11] (TA807853)
ACCN:	NM_001145829 , NP_001139301
Synonyms:	AMCD2B; DA2B; DA2B1; FSSV; fsTnl
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001139301
Locus ID:	7136
Cytogenetics:	11p15.5



[View online »](#)