

## Product datasheet for **TB427498**

### FGF13 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FGF13 (NM_001139500), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC227498
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	FGF13 Mouse Monoclonal Antibody [Clone ID: OTI7E8] (TA807761)
ACCN:	<u><a href="#">NM_001139500</a></u> , <u><a href="#">NP_001132972</a></u>
Synonyms:	DEE90; FGF-13; FGF2; FHF-2; FHF2; LINC00889
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001132972</a></u>
Locus ID:	2258
Cytogenetics:	Xq26.3-q27.1
Protein Families:	Secreted Protein
Protein Pathways:	MAPK signaling pathway, Melanoma, Pathways in cancer, Regulation of actin cytoskeleton



[View online »](#)