

## Product datasheet for **TB427934**

### Periostin (POSTN) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of POSTN (NM_001135935), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC227934
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Periostin (POSTN) Mouse Monoclonal Antibody [Clone ID: OT12B2] (TA804575)
ACCN:	<u><a href="#">NM_001135935</a></u> , <u><a href="#">NP_001129407</a></u>
Synonyms:	OSF-2; OSF2; PDLPOSTN; PN
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001129407</a></u>
Locus ID:	10631
Cytogenetics:	13q13.3
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)