

## Product datasheet for **TB428229**

### Sex Hormone Binding Globulin (SHBG) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SHBG (NM_001146280), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC228229
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Sex Hormone Binding Globulin (SHBG) Mouse Monoclonal Antibody [Clone ID: OT11H10] (TA507187)
ACCN:	<a href="#">NM_001146280</a> , <a href="#">NP_001139752</a>
Synonyms:	ABP; SBP; TEBG
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001139752</a>
Locus ID:	6462
Cytogenetics:	17p13.1
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)