

## Product datasheet for **TB428289**

### Aminomethyltransferase (AMT) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of AMT (NM_001164710), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC228289
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Aminomethyltransferase (AMT) Mouse Monoclonal Antibody [Clone ID: OT16F4] (TA807784)
ACCN:	<a href="#">NM_001164710</a> , <a href="#">NP_001158182</a>
Synonyms:	GCE; GCST; GCVT; NKH
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001158182</a>
Locus ID:	275
Cytogenetics:	3p21.31
Protein Pathways:	Glycine, serine and threonine metabolism, Metabolic pathways, Nitrogen metabolism, One carbon pool by folate



[View online »](#)