

## Product datasheet for **TB428362**

### RRM2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RRM2 (NM_001165931), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC228362
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	RRM2 Mouse Monoclonal Antibody [Clone ID: OT13D3] (TA810467)
ACCN:	<a href="#">NM_001165931</a> , <a href="#">NP_001159403</a>
Synonyms:	C2orf48; R2; RR2; RR2M
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001159403</a>
Locus ID:	6241
Cytogenetics:	2p25.1
Protein Families:	Druggable Genome
Protein Pathways:	Glutathione metabolism, Metabolic pathways, p53 signaling pathway, Purine metabolism, Pyrimidine metabolism



[View online »](#)