

## Product datasheet for **TB429023**

### ELAC2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ELAC2 (NM_173717), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC229023
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_173717</a></u> , <u><a href="#">NP_776065</a></u>
Synonyms:	COXPD17; ELC2; HPC2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_776065</a></u>
Locus ID:	60528
Cytogenetics:	17p12
Protein Families:	Druggable Genome, Stem cell - Pluripotency



[View online »](#)