

Product datasheet for **TB429331**

ZNF135 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ZNF135 (NM_003436), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC229331
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_003436 , NP_003427
Synonyms:	pHZ-17; pT3; ZNF61; ZNF78L1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003427
Locus ID:	7694
Cytogenetics:	19q13.43
Protein Families:	Transcription Factors



[View online »](#)