

Product datasheet for **TB429651**

LOX CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of LOX (NM_001178102), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC229651
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	LOX Mouse Monoclonal Antibody [Clone ID: OTI6B7] (TA500869)
ACCN:	NM_001178102 , NP_001171573
Synonyms:	AAT10
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001171573
Locus ID:	4015
Cytogenetics:	5q23.1
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)