

Product datasheet for **TB429754**

CCRK (CDK20) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CDK20 (NM_001170640), transcript variant 5, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC229754
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CCRK (CDK20) Mouse Monoclonal Antibody [Clone ID: OT11B3] (TA809548)
ACCN:	<u>NM_001170640</u> , <u>NP_001164111</u>
Synonyms:	CCRK; CDCH; P42; PNQALRE
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001164111</u>
Locus ID:	23552
Cytogenetics:	9q22.1
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)