

Product datasheet for **TB429806**

CD244 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CD244 (NM_001166664), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC229806
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CD244 Mouse Monoclonal Antibody [Clone ID: OT11B12] (TA812270)
ACCN:	NM_001166664 , NP_001160136
Synonyms:	2B4; NAIL; NKR2B4; Nmrk; SLAMF4
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001160136
Locus ID:	51744
Cytogenetics:	1q23.3
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Natural killer cell mediated cytotoxicity



[View online »](#)