

## Product datasheet for **TB429981**

### SIGLEC6 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SIGLEC6 (NM_001177548), transcript variant 5, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC229981
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_001177548</a></u> , <u><a href="#">NP_001171019</a></u>
Synonyms:	CD33L; CD33L1; CD33L2; CD327; CDW327; OBBP1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001171019</a></u>
Locus ID:	946
Cytogenetics:	19q13.41
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane



[View online »](#)