

## Product datasheet for **TB429996**

### UDP glucose dehydrogenase (UGDH) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of UGDH (NM_001184701), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC229996
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	UDP glucose dehydrogenase (UGDH) Mouse Monoclonal Antibody [Clone ID: OTI2C11] (TA810531)
ACCN:	<u><a href="#">NM_001184701</a></u> , <u><a href="#">NP_001171630</a></u>
Synonyms:	DEE84; EIEE84; GDH; UDP-GlcDH; UDPGDH; UGD
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001171630</a></u>
Locus ID:	7358
Cytogenetics:	4p14
Protein Pathways:	Amino sugar and nucleotide sugar metabolism, Ascorbate and aldarate metabolism, Metabolic pathways, Pentose and glucuronate interconversions, Starch and sucrose metabolism



[View online »](#)