

## Product datasheet for **TB430242**

### HPSE2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HPSE2 (NM_001166244), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC230242
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_001166244</a></u> , <u><a href="#">NP_001159716</a></u>
Synonyms:	HPA2; HPR2; UFS; UFS1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001159716</a></u>
Locus ID:	60495
Cytogenetics:	10q24.2
Protein Pathways:	Glycosaminoglycan degradation, Metabolic pathways



[View online »](#)