

Product datasheet for **TB430385**

XPNPEP1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of XPNPEP1 (NM_001167604), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC230385
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	XPNPEP1 Mouse Monoclonal Antibody [Clone ID: OTI1E3] (TA501383)
ACCN:	NM_001167604 , NP_001161076
Synonyms:	APP1; SAMP; XPNPEP; XPNPEPL; XPNPEPL1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001161076
Locus ID:	7511
Cytogenetics:	10q25.1
Protein Families:	Druggable Genome, Protease



[View online »](#)