

## Product datasheet for **TB430457**

### BCAR1 CytoSection

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | CytoSections  |
| Description:                          | Transient overexpression of BCAR1 (NM_001170720), transcript variant 8, in HEK293T cells, paraffin embedded controls for ICC/IHC staining   |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC230457  |
| Tag:                                  | C-MYC/DDK   |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)  |
| Target Detection Antibodies:          | BCAR1 Mouse Monoclonal Antibody [Clone ID: OT11A9] (TA502613)   |
| ACCN:                                 | <u><a href="#">NM_001170720</a></u> , <u><a href="#">NP_001164191</a></u>   |
| Synonyms:                             | CAS; CAS1; CASS1; CRKAS; P130Cas  |
| Storage:                              | Room Temperature, or 2-8°C for long term storage  |
| Stability:                            | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| RefSeq:                               | <u><a href="#">NP_001164191</a></u>   |
| Locus ID:                             | 9564  |
| Cytogenetics:                         | 16q23.1   |
| Protein Families:                     | Druggable Genome  |
| Protein Pathways:                     | Chemokine signaling pathway, Focal adhesion, Leukocyte transendothelial migration, Regulation of actin cytoskeleton   |



[View online »](#)