

## Product datasheet for **TB430550**

### BCAR1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of BCAR1 (NM_001170719), transcript variant 7, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC230550
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	BCAR1 Mouse Monoclonal Antibody [Clone ID: OT11A9] (TA502613)
ACCN:	<a href="#">NM_001170719</a> , <a href="#">NP_001164190</a>
Synonyms:	CAS; CAS1; CASS1; CRKAS; P130Cas
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001164190</a>
Locus ID:	9564
Cytogenetics:	16q23.1
Protein Families:	Druggable Genome
Protein Pathways:	Chemokine signaling pathway, Focal adhesion, Leukocyte transendothelial migration, Regulation of actin cytoskeleton



[View online »](#)