

## Product datasheet for **TB430611**

### GPR64 (ADGRG2) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ADGRG2 (NM_001184833), transcript variant 5, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC230611
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001184833</u> , <u>NP_001171762</u>
Synonyms:	CBAVDX; EDDM6; GPR64; HE6; TM7LN2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001171762</u>
Locus ID:	10149
Cytogenetics:	Xp22.13
Protein Families:	Druggable Genome, GPCR, Transmembrane



[View online »](#)