

Product datasheet for **TB430941**

TRAIL (TNFSF10) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TNFSF10 (NM_001190942), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC230941
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TRAIL (TNFSF10) Mouse Monoclonal Antibody [Clone ID: OT112E12] (TA811771)
ACCN:	NM_001190942 , NP_001177871
Synonyms:	Apo-2L; APO2L; CD253; TL2; TNLG6A; TRAIL
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001177871
Locus ID:	8743
Cytogenetics:	3q26.31
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Apoptosis, Cytokine-cytokine receptor interaction, Natural killer cell mediated cytotoxicity



[View online »](#)