

## Product datasheet for **TB431228**

### Dystrobrevin alpha (DTNA) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DTNA (NM_001198943), transcript variant 15, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC231228
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Dystrobrevin alpha (DTNA) Mouse Monoclonal Antibody [Clone ID: OT11B2] (TA502184)
ACCN:	<a href="#">NM_001198943</a> , <a href="#">NP_001185872</a>
Synonyms:	D18S892E; DRP3; DTN; DTN-A; LVNC1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001185872</a>
Locus ID:	1837
Cytogenetics:	18q12.1
Protein Families:	Druggable Genome



[View online »](#)