

Product datasheet for **TB431368**

CRMP2 (DPYSL2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DPYSL2 (NM_001197293), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC231368
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CRMP2 (DPYSL2) Mouse Monoclonal Antibody [Clone ID: OTI5D7] (TA811601)
ACCN:	NM_001197293 , NP_001184222
Synonyms:	CRMP-2; CRMP2; DHPRP2; DRP-2; DRP2; N2A3; ULIP-2; ULIP2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001184222
Locus ID:	1808
Cytogenetics:	8p21.2
Protein Families:	Druggable Genome
Protein Pathways:	Axon guidance



[View online »](#)