

## Product datasheet for **TB431372**

### LDL Receptor (LDLR) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of LDLR (NM_001195803), transcript variant 6, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC231372
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	LDL Receptor (LDLR) Mouse Monoclonal Antibody [Clone ID: OT17C5] (TA812510)
ACCN:	<a href="#">NM_001195803</a> , <a href="#">NP_001182732</a>
Synonyms:	FH; FHC; FHCL1; LDLCQ2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001182732</a>
Locus ID:	3949
Cytogenetics:	19p13.2
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Endocytosis



[View online »](#)