

Product datasheet for **TB431901**

Nkx3.1 (NKX3-1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NKX3-1 (NM_001256339), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC231901
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Nkx3.1 (NKX3-1) Mouse Monoclonal Antibody [Clone ID: OTI6E7] (TA805425)
ACCN:	NM_001256339 , NP_001243268
Synonyms:	BAPX2; NKX3; NKX3.1; NKX3A
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001243268
Locus ID:	4824
Cytogenetics:	8p21.2
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Pathways in cancer, Prostate cancer



[View online »](#)