

## Product datasheet for **TB432050**

### ETV7 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ETV7 (NM_001207041), transcript variant 8, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC232050
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ETV7 Mouse Monoclonal Antibody [Clone ID: OTI3B2] (TA802936)
ACCN:	<a href="#">NM_001207041</a> , <a href="#">NP_001193970</a>
Synonyms:	TEL-2; TEL2; TELB
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001193970</a>
Locus ID:	51513
Cytogenetics:	6p21.31
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Dorso-ventral axis formation



[View online »](#)