

Product datasheet for **TB432178**

INPP5F CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of INPP5F (NM_001243195), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC232178
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	INPP5F Mouse Monoclonal Antibody [Clone ID: OT11D4] (TA809925)
ACCN:	NM_001243195 , NP_001230124
Synonyms:	hSAC2; MSTP007; MSTPO47; SAC2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001230124
Locus ID:	22876
Cytogenetics:	10q26.11
Protein Families:	Druggable Genome



[View online »](#)