

Product datasheet for **TB432250**

CD63 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CD63 (NM_001257389), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC232250
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CD63 Mouse Monoclonal Antibody [Clone ID: OTI2G6] (TA802751)
ACCN:	NM_001257389 , NP_001244318
Synonyms:	LAMP-3; ME491; MLA1; OMA81H; TSPAN30
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001244318
Locus ID:	967
Cytogenetics:	12q13.2
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Lysosome



[View online »](#)