

## Product datasheet for **TB432507**

### BTN3A2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of BTN3A2 (NM_001197246), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC232507
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	BTN3A2 Mouse Monoclonal Antibody [Clone ID: OT11A6] (TA500730)
ACCN:	<a href="#">NM_001197246</a> , <a href="#">NP_001184175</a>
Synonyms:	BT3.2; BTF4; BTN3.2; CD277
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001184175</a>
Locus ID:	11118
Cytogenetics:	6p22.2
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)