

Product datasheet for **TB432525**

Progesterone Receptor (PGR) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PGR (NM_001271162), transcript variant 7, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC232525
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Progesterone Receptor (PGR) Mouse Monoclonal Antibody [Clone ID: OT111E8] (TA802704)
ACCN:	NM_001271162 , NP_001258091
Synonyms:	NR3C3; PR
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001258091
Locus ID:	5241
Cytogenetics:	11q22.1
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors
Protein Pathways:	Oocyte meiosis, Progesterone-mediated oocyte maturation



[View online »](#)