

## Product datasheet for **TB432632**

### Bif (SH3GLB1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SH3GLB1 (NM_001206652), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC232632
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001206652</u> , <u>NP_001193581</u>
Synonyms:	Bif-1; CGI-61; dj612B15.2; PPP1R70
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001193581</u>
Locus ID:	51100
Cytogenetics:	1p22.3
Protein Pathways:	Endocytosis



[View online »](#)