

Product datasheet for **TB432633**

GDAP1L1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GDAP1L1 (NM_001256737), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC232633
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	GDAP1L1 Mouse Monoclonal Antibody [Clone ID: OT11G5] (TA503153)
ACCN:	NM_001256737 , NP_001243666
Synonyms:	dj881L22.1; dj995J12.1.1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001243666
Locus ID:	78997
Cytogenetics:	20q13.12
Protein Families:	Transmembrane



[View online »](#)