

Product datasheet for **TB432754**

TAF1 48 (TAF1A) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TAF1A (NM_001201536), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC232754
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TAF1 48 (TAF1A) Mouse Monoclonal Antibody [Clone ID: OTI7A11] (TA811372)
ACCN:	NM_001201536 , NP_001188465
Synonyms:	MGC:17061; RAFI48; SL1; TAFI48
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001188465
Locus ID:	9015
Cytogenetics:	1q41
Protein Families:	Transcription Factors



[View online »](#)