

Product datasheet for **TB433006**

Acetyl CoA synthetase (ACSS2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ACSS2 (NM_001242393), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC233006
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Acetyl CoA synthetase (ACSS2) Mouse Monoclonal Antibody [Clone ID: OTI2C6] (TA503610)
ACCN:	NM_001242393 , NP_001229322
Synonyms:	ACAS2; ACECS; AceCS1; ACS; ACSA; dj1161H23.1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001229322
Locus ID:	55902
Cytogenetics:	20q11.22
Protein Pathways:	Glycolysis / Gluconeogenesis, Metabolic pathways, Propanoate metabolism, Pyruvate metabolism



[View online »](#)