

Product datasheet for **TB433029**

Transketolase (TKT) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TKT (NM_001258028), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC233029
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Transketolase (TKT) Mouse Monoclonal Antibody [Clone ID: OTI5H3] (TA500892)
ACCN:	NM_001258028 , NP_001244957
Synonyms:	HEL-S-48; HEL107; SDDHD; TK; TKT1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001244957
Locus ID:	7086
Cytogenetics:	3p21.1
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Pentose phosphate pathway



[View online »](#)