

Product datasheet for **TB433304**

Prostaglandin dehydrogenase 1 (HPGD) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HPGD (NM_001256307), transcript variant 6, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC233304
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Prostaglandin dehydrogenase 1 (HPGD) Mouse Monoclonal Antibody [Clone ID: OTI2C10] (TA502561)
ACCN:	NM_001256307 , NP_001243236
Synonyms:	15-PGDH; PGDH; PGDH1; PHOAR1; SDR36C1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001243236
Locus ID:	3248
Cytogenetics:	4q34.1
Protein Families:	Druggable Genome



[View online »](#)