

Product datasheet for **TB433516**

EMA (MUC1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MUC1 (NM_001204292), transcript variant 16, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC233516
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	EMA (MUC1) Mouse Monoclonal Antibody [Clone ID: OT12F6] (TA800790)
ACCN:	<u>NM_001204292</u> , <u>NP_001191221</u>
Synonyms:	ADMCKD; ADMCKD1; ADTKD2; CA 15-3; CD227; EMA; H23AG; KL-6; MAM6; MCD; MCKD; MCKD1; MUC-1; MUC-1/SEC; MUC-1/X; MUC1/ZD; PEM; PEMT; PUM
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001191221</u>
Locus ID:	4582
Cytogenetics:	1q22
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane



[View online »](#)