

Product datasheet for **TB433700**

Aurora B (AURKB) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of AURKB (NM_001256834), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC233700
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001256834 , NP_001243763
Synonyms:	AIK2; AIM-1; AIM1; ARK-2; ARK2; AurB; aurkb-sv1; aurkb-sv2; IPL1; PPP1R48; STK-1; STK5; STK12
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001243763
Locus ID:	9212
Cytogenetics:	17p13.1
Protein Families:	Druggable Genome, Protein Kinase, Stem cell - Pluripotency



[View online »](#)