

Product datasheet for **TB433709**

GIRK1 (KCNJ3) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KCNJ3 (NM_001260509), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC233709
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	GIRK1 (KCNJ3) Mouse Monoclonal Antibody [Clone ID: OTI3E11] (TA504152)
ACCN:	NM_001260509 , NP_001247438
Synonyms:	GIRK1; KGA; KIR3.1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001247438
Locus ID:	3760
Cytogenetics:	2q24.1
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane



[View online »](#)